

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.82	
Mean	176.55	Mean	534.65
SD	3.29	SD	9.73
CV	1.86	CV	1.82
Day-Run	Value	Day-Run	Value
1 - 1	178.2	1 - 1	551
2 - 1	179	2 - 1	540
7 - 1	172.8	7 - 1	528.8
8 - 1	180	8 - 1	545
9 - 1	176.1	9 - 1	529.1
11 - 1	175	11 - 1	514
12 - 1	181.9	12 - 1	530
14 - 1	179	14 - 1	530
15 - 1	180.4	15 - 1	541
16 - 1	176	16 - 1	517.9
17 - 1	178	17 - 1	538
18 - 1	176	18 - 1	525.3
19 - 1	180	19 - 1	535
21 - 1	178.5	21 - 1	539.5
22 - 1	178.5	22 - 1	545
23 - 1	178	23 - 1	552
24 - 1	172.8	24 - 1	527.5
25 - 1	174	25 - 1	542
26 - 1	170.3	26 - 1	532.5
28 - 1	171	28 - 1	538.3
29 - 1	177	29 - 1	526
30 - 1	171.5	30 - 1	534.4

Westgard	تفسیر	Ikan Phosphatase(ALP)	آزمایش	اتوآنالایزر BS-800	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	176.13	SD	2.66	Exp. Date	1403/1/17	Lot. No	035

