

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.37	
Mean	29.01	Mean	136.68
SD	1.16	SD	3.23
CV	3.98	CV	2.37
Day-Run	Value	Day-Run	Value
1 - 1	29	1 - 1	136.2
3 - 1	28.7	3 - 1	140.6
4 - 1	28.4	4 - 1	137.5
5 - 1	30.2	5 - 1	132.8
6 - 1	28.8	6 - 1	139.5
10 - 1	30.2	10 - 1	132.7
11 - 1	30.4	11 - 1	133
13 - 1	29	13 - 1	136.2
19 - 1	30.3	19 - 1	143
20 - 1	29	20 - 1	132.6
21 - 1	26.8	21 - 1	134.9
24 - 1	30.3	24 - 1	137
25 - 1	28.6	25 - 1	136.8
26 - 1	30.5	26 - 1	138.4
27 - 1	29.2	27 - 1	141
28 - 1	27	28 - 1	140
29 - 1	27.7	29 - 1	133
31 - 1	28	31 - 1	135

Westgard	تفسیر	SGPT (ALT)	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Normal Control	biorex	Mean	29.07	SD	1.02	Exp. Date	1403/1/17	Lot. No	035

