

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 1.57	
Mean	27.9	Mean	135.75
SD	0.81	SD	2.13
CV	2.9	CV	1.57
Day-Run	Value	Day-Run	Value
1 - 1	28.2	1 - 1	136.9
2 - 1	27.9	2 - 1	138.1
7 - 1	28.4	7 - 1	139.1
8 - 1	27	8 - 1	136.1
9 - 1	28.7	9 - 1	137
11 - 1	29.1	11 - 1	134.8
12 - 1	27	12 - 1	132.9
14 - 1	28	14 - 1	132.7
15 - 1	29.4	15 - 1	137
16 - 1	27	16 - 1	135.6
17 - 1	28.6	17 - 1	131.7
18 - 1	26.6	18 - 1	133.7
19 - 1	28.2	19 - 1	135
21 - 1	27.8	21 - 1	136.3
22 - 1	28.7	22 - 1	137
23 - 1	28.7	23 - 1	134.5
24 - 1	27.4	24 - 1	137.3
25 - 1	28	25 - 1	140
26 - 1	28	26 - 1	135.6
28 - 1	26.5	28 - 1	134.8
29 - 1	27	29 - 1	133
30 - 1	27.7	30 - 1	137.3

Westgard	تفسیر	SGPT ( ALT )	آزمایش	اتوآنالایزر BS-800	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	28.01	SD	1.3	Exp. Date	1403/1/17	Lot. No	035

