

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.38	
Mean	38.65	Mean	143.37
SD	0.99	SD	1.97
CV	2.57	CV	1.38
Day-Run	Value	Day-Run	Value
1 - 1	38.8	1 - 1	142.5
3 - 1	40.3	3 - 1	145
4 - 1	39	4 - 1	142.9
5 - 1	37	5 - 1	144.5
6 - 1	39.4	6 - 1	144.9
10 - 1	38	10 - 1	144.7
11 - 1	38.4	11 - 1	140.5
13 - 1	37.7	13 - 1	143.4
19 - 1	39.7	19 - 1	146.8
20 - 1	39.2	20 - 1	142.4
21 - 1	39.5	21 - 1	142
24 - 1	39.9	24 - 1	138.8
25 - 1	38	25 - 1	144.6
26 - 1	39.6	26 - 1	143
27 - 1	38.6	27 - 1	143.7
28 - 1	38	28 - 1	142.5
29 - 1	37.6	29 - 1	146.4
31 - 1	37	31 - 1	142

Westgard	تفسیر	SGOT (AST)	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Normal Control	biorex	Mean	38.04	SD	1.2	Exp. Date	1403/1/17	Lot. No	035

