

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	38.81
SD	0.91
CV	2.34
Day-Run	Value
1 - 1	37.5
2 - 1	37.9
7 - 1	37.7
8 - 1	38.5
9 - 1	39.1
11 - 1	39.4
12 - 1	38.7
14 - 1	39.8
15 - 1	39.9
16 - 1	38.8
17 - 1	40
18 - 1	38.3
19 - 1	38.3
21 - 1	38.1
22 - 1	38.7
23 - 1	39.5
24 - 1	37.4
25 - 1	38.4
26 - 1	39.7
28 - 1	40.3
29 - 1	37.8
30 - 1	40.1

High - CV : 1.55	
Mean	144.1
SD	2.23
CV	1.55
Day-Run	Value
1 - 1	144.5
2 - 1	146.2
7 - 1	148.1
8 - 1	146.1
9 - 1	147.6
11 - 1	146
12 - 1	146.6
14 - 1	144.7
15 - 1	145.7
16 - 1	142
17 - 1	143.4
18 - 1	144
19 - 1	144.8
21 - 1	141
22 - 1	142
23 - 1	143.8
24 - 1	141.1
25 - 1	141.8
26 - 1	141.3
28 - 1	141
29 - 1	145.6
30 - 1	142.8

Westgard	تفسیر	SGOT ( AST )	آزمایش	اتوآنالایزر BS-800	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	38.65	SD	0.99	Exp. Date	1403/1/17	Lot. No	035

