

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

| Normal | | High - CV : 2.15 | |
|---------|-------|------------------|-------|
| Mean | 4.28 | Mean | 4.92 |
| SD | 0.11 | SD | 0.11 |
| CV | 2.48 | CV | 2.15 |
| Day-Run | Value | Day-Run | Value |
| 1 - 1 | 4.4 | 1 - 1 | 4.9 |
| 3 - 1 | 4.1 | 3 - 1 | 5.1 |
| 4 - 1 | 4.3 | 4 - 1 | 5 |
| 5 - 1 | 4.4 | 5 - 1 | 5 |
| 6 - 1 | 4.3 | 6 - 1 | 5 |
| 10 - 1 | 4.3 | 10 - 1 | 5 |
| 11 - 1 | 4.2 | 11 - 1 | 4.8 |
| 13 - 1 | 4.4 | 13 - 1 | 5 |
| 19 - 1 | 4.3 | 19 - 1 | 5 |
| 20 - 1 | 4.3 | 20 - 1 | 4.9 |
| 21 - 1 | 4.3 | 21 - 1 | 5 |
| 24 - 1 | 4.4 | 24 - 1 | 4.8 |
| 25 - 1 | 4.3 | 25 - 1 | 4.9 |
| 26 - 1 | 4.4 | 26 - 1 | 5 |
| 27 - 1 | 4.1 | 27 - 1 | 4.8 |
| 28 - 1 | 4.1 | 28 - 1 | 4.7 |
| 29 - 1 | 4.2 | 29 - 1 | 4.9 |
| 31 - 1 | 4.2 | 31 - 1 | 4.8 |

| | | | | | | | | | |
|----------------|--------|---------|--------|--------------------|--------|-----------|-----------|---------|-----|
| Westgard | تفسیر | Albomin | آزمایش | BS-800 اتوآنالایزر | دستگاه | 4 | ماه | 1402 | سال |
| Normal Control | biorex | Mean | 4.26 | SD | 0.15 | Exp. Date | 1403/1/17 | Lot. No | 035 |

