

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

| Normal | | High - CV : 2.13 | |
|---------|--------|------------------|--------|
| Mean | 176.13 | Mean | 538.65 |
| SD | 2.66 | SD | 11.48 |
| CV | 1.51 | CV | 2.13 |
| Day-Run | Value | Day-Run | Value |
| 1 - 1 | 177 | 1 - 1 | 524.5 |
| 3 - 1 | 179.8 | 3 - 1 | 538.5 |
| 4 - 1 | 171 | 4 - 1 | 535.8 |
| 5 - 1 | 176 | 5 - 1 | 552.1 |
| 6 - 1 | 176.8 | 6 - 1 | 548.1 |
| 10 - 1 | 176.8 | 10 - 1 | 535.7 |
| 11 - 1 | 179 | 11 - 1 | 541 |
| 13 - 1 | 170 | 13 - 1 | 547.3 |
| 19 - 1 | 173.1 | 19 - 1 | 543.7 |
| 20 - 1 | 174.3 | 20 - 1 | 528.7 |
| 21 - 1 | 176 | 21 - 1 | 526.2 |
| 24 - 1 | 178 | 24 - 1 | 522.1 |
| 25 - 1 | 180 | 25 - 1 | 533 |
| 26 - 1 | 177 | 26 - 1 | 530 |
| 27 - 1 | 177 | 27 - 1 | 544 |
| 28 - 1 | 176.6 | 28 - 1 | 537 |
| 29 - 1 | 175.6 | 29 - 1 | 570 |
| 31 - 1 | 176.4 | 31 - 1 | 538 |

| | | | | | | | | | |
|----------------|--------|-----------------------|--------|--------------------|--------|-----------|-----------|---------|-----|
| Westgard | تفسیر | Ikan Phosphatase(ALP) | آزمایش | اتوآنالایزر BS-800 | دستگاه | 4 | ماه | 1402 | سال |
| Normal Control | biorex | Mean | 175.61 | SD | 3.87 | Exp. Date | 1403/1/17 | Lot. No | 035 |

