

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	79	
3 - 1	80.7	
4 - 1	82.6	
5 - 1	80.1	
6 - 1	80.1	
10 - 1	80.1	
11 - 1	83	
13 - 1	81.1	
19 - 1	81	
20 - 1	79.8	
21 - 1	77	
21 - 1	77	
24 - 1	83	
25 - 1	77.5	
26 - 1	82	
27 - 1	79.4	
28 - 1	82.2	
29 - 1	80	
31 - 1	79	

Westgard	تفسیر	Amylase	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	79.25	SD	2.2	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	80.24	SD	1.85	CV	2.31			

