

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	80.5	
2 - 1	78	
7 - 1	76.4	1:2S
8 - 1	80.5	
9 - 1	79	
11 - 1	82	1:2S
12 - 1	76.5	
14 - 1	83	1:2S
15 - 1	77	
16 - 1	81.4	
17 - 1	83.1	1:2S
18 - 1	77.8	
19 - 1	82	1:2S
21 - 1	77	
22 - 1	80.4	
23 - 1	76.8	
24 - 1	81.4	
25 - 1	81.3	
26 - 1	80	
28 - 1	78.4	
29 - 1	81.7	
30 - 1	81	

Westgard	تفسیر	Amylase	آزمایش	BS-800 اتوآنالایزر	دستگاه	5	ماه	1402	سال
Control	biorex	Mean	78.24	SD	1.85	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	79.78	SD	2.2	CV	2.76			

