

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	20.58	
2 - 1	21.3	
3 - 1	20	
4 - 1	20.65	
7 - 1	21	
8 - 1	20.61	
9 - 1	21	
10 - 1	21.25	
11 - 1	20.18	
12 - 1	20.54	
14 - 1	20.95	
15 - 1	20.41	
16 - 1	20.32	
17 - 1	20.52	
18 - 1	21.1	
19 - 1	21.06	Cusum 3S
21 - 1	20.95	
22 - 1	20.22	
23 - 1	20.12	
24 - 1	20.59	
25 - 1	20.65	
26 - 1	20.71	
28 - 1	20.5	
29 - 1	21.15	
30 - 1	21.2	
31 - 1	20.53	

Westgard	تفسیر	BHCG	آزمایش	2000 ایمولایت XPI-	دستگاه	5	ماه	1402	سال
Control	biorex	Mean	20.56	SD	0.72	Exp. Date	1402/10/30	Lot. No	36-02-028
Next Month	Mean	20.7	SD	0.37	CV	1.79			

