

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	1
SD	0.06
CV	5.74
Day-Run	Value
1 - 1	0.96
2 - 1	1.03
2 - 1	1.03
7 - 1	0.98
8 - 1	0.91
9 - 1	1.05
11 - 1	1
12 - 1	1.07
14 - 1	0.91
15 - 1	1
16 - 1	0.97
17 - 1	1.02
18 - 1	0.95
19 - 1	0.92
21 - 1	1.06
22 - 1	0.98
23 - 1	0.92
24 - 1	0.95
24 - 1	0.95
25 - 1	1.03
26 - 1	1.1
28 - 1	1.04
29 - 1	1.1
30 - 1	1

High - CV : 2.9	
Mean	4.26
SD	0.12
CV	2.9
Day-Run	Value
1 - 1	4.25
2 - 1	4.25
7 - 1	4.3
8 - 1	4.46
9 - 1	4.11
11 - 1	4.3
11 - 2	4.3
12 - 1	4.11
14 - 1	4.5
15 - 1	4.11
16 - 1	4.11
17 - 1	4.42
18 - 1	4.05
19 - 1	4.4
21 - 1	4.17
22 - 1	4.37
23 - 1	4.15
24 - 1	4.25
25 - 1	4.34
25 - 2	4.34
26 - 1	4.19
28 - 1	4.36
29 - 1	4.18
30 - 1	4.32

Westgard	تفسیر	Bili T	آزمایش	اتوآنالایزر BS-800	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	0.98	SD	0.06	Exp. Date	1403/1/17	Lot. No	035

