

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	80	
3 - 1	84.2	
4 - 1	87	
6 - 1	82.4	
7 - 1	84.5	
10 - 1	86	
11 - 1	88	
12 - 1	82.6	
13 - 1	83.3	
14 - 1	83.6	
15 - 2	87.3	
17 - 1	85.2	
18 - 1	83.9	
19 - 1	85	
20 - 1	80.4	
21 - 1	86.6	
22 - 1	87	
24 - 1	86.8	
25 - 1	89.2	
26 - 1	87.8	4:1S
27 - 1	85.2	
28 - 1	84	
29 - 1	87.5	
31 - 1	86.3	

Westgard	تفسیر	C3	آزمایش	BS-480 اتوآنالایزر	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	84.39	SD	2.35	Exp. Date	1402/5/1	Lot. No	641-01-030/2
Next Month	Mean	85.16	SD	2.36	CV	2.78			

