

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	48.68	
2 - 1	49.95	
3 - 1	51.2	
4 - 1	50.01	
7 - 1	50.89	
8 - 1	50.2	
9 - 1	51.5	
10 - 1	49.7	
11 - 1	51.8	
12 - 1	51.7	
14 - 1	49.8	
15 - 1	47.6	
16 - 1	48.5	
17 - 1	50.1	
18 - 1	51.58	
19 - 1	51.03	
21 - 1	50.86	Cusum 3S
22 - 1	46.5	
23 - 1	47.9	
24 - 1	47.2	
25 - 1	48.1	
26 - 1	51.8	
28 - 1	48.9	
29 - 1	51.1	
30 - 1	51.8	
31 - 1	50.7	

Westgard	تفسیر	CA 125	آزمایش	2000 ایمولایت XPI-	دستگاه	5	ماه	1402	سال
Control	Seronorm	Mean	49.88	SD	2.29	Exp. Date	1402/10/30	Lot. No	2011909
Next Month	Mean	49.97	SD	1.58	CV	3.16			

