

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	18	
3 - 1	17.5	
6 - 1	17	
8 - 1	17.2	
10 - 1	18.5	
13 - 1	19	1:2S
16 - 1	17.9	
17 - 1	19.1	1:2S
20 - 1	17.5	
22 - 1	18	
27 - 1	17.8	
29 - 1	18.5	
31 - 1	17.8	

Westgard	تفسیر	CEA	آزمایش	SIEMENS Centaur CP	دستگاه	3	ماه	1402	سال
Control	biorex	Mean	17.67	SD	0.65	Exp. Date	1402/4/19	Lot. No	
Next Month	Mean	17.98	SD	0.64	CV	3.57			

