

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	126.68
SD	1.92
CV	1.51
Day-Run	Value
1 - 1	126
3 - 1	126
3 - 1	126
4 - 1	129
5 - 1	129
6 - 1	127
10 - 1	128
11 - 1	129
13 - 1	129
19 - 1	129
20 - 1	127
21 - 1	123
24 - 1	126
25 - 1	124
26 - 1	127
27 - 1	125
28 - 1	128
29 - 1	125
31 - 1	124

High - CV : 2.01	
Mean	212.94
SD	4.28
CV	2.01
Day-Run	Value
1 - 1	210
3 - 1	215
4 - 1	220
5 - 1	211
6 - 1	213
10 - 1	212
11 - 1	206
13 - 1	210
19 - 1	216
20 - 1	210
21 - 1	205
24 - 1	216
25 - 1	210
26 - 1	219
27 - 1	213
28 - 1	220
29 - 1	213
31 - 1	214

Westgard	تفسیر	CPK	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Normal Control	biorex	Mean	126.9	SD	1.8	Exp. Date	1403/1/17	Lot. No	035

