

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.48	
Mean	126.9	Mean	214.25
SD	1.8	SD	3.18
CV	1.42	CV	1.48
Day-Run	Value	Day-Run	Value
1 - 1	126	1 - 1	213
2 - 1	125	2 - 1	216
3 - 1	124	3 - 1	211
6 - 1	130	6 - 1	219
8 - 1	128	8 - 1	219
9 - 1	129	9 - 1	218
10 - 1	128	10 - 1	213
11 - 1	130	11 - 1	212
17 - 1	129	17 - 1	217
20 - 1	127	20 - 1	215
21 - 1	128	21 - 1	211
22 - 1	125	22 - 1	215
23 - 1	125	23 - 1	217
24 - 1	126	24 - 1	209
25 - 1	127	25 - 1	210
27 - 1	124	27 - 1	217
28 - 1	127	28 - 1	210
29 - 1	126	29 - 1	214
30 - 1	127	30 - 1	212
31 - 1	127	31 - 1	217

Westgard	تفسیر	CPK	آزمایش	اتوآنالایزر BS-800	دستگاه	3	ماه	1402	سال
Normal Control	biorex	Mean	127.07	SD	1.5	Exp. Date	1403/1/17	Lot. No	035

