

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.31	
Mean	9.46	Mean	14.09
SD	0.2	SD	0.33
CV	2.12	CV	2.31
Day-Run	Value	Day-Run	Value
1 - 1	9.8	1 - 1	14
3 - 1	9.7	3 - 1	14.7
4 - 1	9.5	4 - 1	14.5
5 - 1	9.7	5 - 1	14.3
6 - 1	9.5	6 - 1	14.5
10 - 1	9.4	10 - 1	14.2
11 - 1	9.2	11 - 1	14
13 - 1	9.5	13 - 1	14.6
19 - 1	9.6	19 - 1	13.7
20 - 1	9.3	20 - 1	14.2
21 - 1	9.2	21 - 1	13.7
24 - 1	9.6	24 - 1	13.8
25 - 1	9.3	25 - 1	14.1
26 - 1	9.5	26 - 1	13.9
27 - 1	9.2	27 - 1	14
28 - 1	9.2	28 - 1	13.8
29 - 1	9.3	29 - 1	14.1
31 - 1	9.7	31 - 1	13.6

Westgard	تفسیر	Calcium	آزمایش	اتوانالایزر BS-800	دستگاه	4	ماه	1402	سال
Normal Control	biorex	Mean	9.49	SD	0.17	Exp. Date	1403/1/17	Lot. No	035

