

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.67	
Mean	9.41	Mean	13.91
SD	0.15	SD	0.23
CV	1.58	CV	1.67
Day-Run	Value	Day-Run	Value
1 - 1	9.2	1 - 1	13.8
2 - 1	9.3	2 - 1	13.7
7 - 1	9.6	7 - 1	14.2
8 - 1	9.5	8 - 1	13.9
9 - 1	9.2	9 - 1	14.1
11 - 1	9.3	11 - 1	13.9
12 - 1	9.6	12 - 1	13.9
14 - 1	9.4	14 - 1	14.3
15 - 1	9.4	15 - 1	13.9
16 - 1	9.5	16 - 1	14
17 - 1	9.3	17 - 1	13.8
18 - 1	9.5	18 - 1	14
19 - 1	9.3	19 - 1	13.8
21 - 1	9.2	21 - 1	14.1
22 - 1	9.4	22 - 1	14.3
23 - 1	9.4	23 - 1	13.7
24 - 1	9.3	24 - 1	13.5
25 - 1	9.6	25 - 1	13.8
26 - 1	9.7	26 - 1	13.5
28 - 1	9.5	28 - 1	14.2
29 - 1	9.3	29 - 1	13.6
30 - 1	9.6	30 - 1	14

Westgard	تفسیر	Calcium	آزمایش	BS-800 اتوآنالایزر	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	9.4	SD	0.25	Exp. Date	1403/1/17	Lot. No	035

