

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 0.93	
Mean	173.02	Mean	262.95
SD	2.6	SD	2.44
CV	1.5	CV	0.93
Day-Run	Value	Day-Run	Value
1 - 1	175	1 - 1	263
2 - 1	170.5	2 - 1	260
3 - 1	174	3 - 1	263
6 - 1	169.8	6 - 1	259
8 - 1	175	8 - 1	263
9 - 1	177	9 - 1	263
10 - 1	176	10 - 1	266
11 - 1	171	11 - 1	260
17 - 1	170	17 - 1	262
20 - 1	174	20 - 1	263
21 - 1	172	21 - 1	260
22 - 1	173	22 - 1	266
23 - 1	173	23 - 1	260
24 - 1	173	24 - 1	266
25 - 1	175	25 - 1	266
27 - 1	170	27 - 1	260
28 - 1	170	28 - 1	264
29 - 1	172	29 - 1	264
30 - 1	171	30 - 1	266
31 - 1	179	31 - 1	265

Westgard	تفسیر	Cholesterol	آزمایش	BS-800 اتوآنالایزر	دستگاه	3	ماه	1402	سال
Normal Control	biorex	Mean	172.89	SD	2.89	Exp. Date	1403/1/17	Lot. No	035

