

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 0.94	
Mean	176	Mean	264.33
SD	2.74	SD	2.5
CV	1.56	CV	0.94
Day-Run	Value	Day-Run	Value
1 - 1	173	1 - 1	267
3 - 1	179	3 - 1	267
4 - 1	176	4 - 1	263
5 - 1	176	5 - 1	265
6 - 1	173	6 - 1	265
10 - 1	178	10 - 1	264
11 - 1	178	11 - 1	265
13 - 1	178	13 - 1	266
19 - 1	174	19 - 1	259
20 - 1	175	20 - 1	265
21 - 1	178	21 - 1	265
24 - 1	171	24 - 1	267
25 - 1	177	25 - 1	263
26 - 1	178	26 - 1	264
27 - 1	177	27 - 1	260
28 - 1	170	28 - 1	266
29 - 1	179	29 - 1	260
31 - 1	178	31 - 1	267

Westgard	تفسیر	Cholesterol	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Normal Control	biorex	Mean	175.02	SD	2.6	Exp. Date	1403/1/17	Lot. No	035

