

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	120	
2 - 1	123.7	
3 - 1	126	
4 - 1	121	
7 - 1	125.5	
8 - 1	125.7	
9 - 1	122	
12 - 1	126	
14 - 1	128	
15 - 1	121.4	
16 - 1	127	
17 - 1	125	
18 - 1	128	Cusum 3S
19 - 1	130	1:2S Cusum 3S
21 - 1	120	
22 - 1	127	
23 - 1	122	
24 - 1	126	
25 - 1	129	1:2S
26 - 1	127	
28 - 1	122	
29 - 1	124	
30 - 1	126	

Westgard	تفسیر	Copper ( Cu )	آزمایش	BS-800 اتوآنالایزر	دستگاه	5	ماه	1402	سال
Control	biorex	Mean	124.04	SD	2.08	Exp. Date	1402/9/20	Lot. No	hsn.35
Next Month	Mean	124.88	SD	2.89	CV	2.32			

