

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.07	
Mean	2.52	Mean	3.61
SD	0.06	SD	0.07
CV	2.28	CV	2.07
Day-Run	Value	Day-Run	Value
1 - 1	2.48	1 - 1	3.56
2 - 1	2.48	2 - 1	3.61
3 - 1	2.4	3 - 1	3.7
6 - 1	2.5	6 - 1	3.7
8 - 1	2.55	8 - 1	3.6
9 - 1	2.45	9 - 1	3.6
10 - 1	2.5	10 - 1	3.47
11 - 1	2.5	11 - 1	3.45
17 - 1	2.54	17 - 1	3.7
21 - 1	2.6	21 - 1	3.6
22 - 1	2.61	22 - 1	3.6
23 - 1	2.5	23 - 1	3.7
24 - 1	2.54	24 - 1	3.54
25 - 1	2.58	25 - 1	3.6
27 - 1	2.57	27 - 1	3.55
28 - 1	2.54	28 - 1	3.65
29 - 1	2.45	29 - 1	3.63
30 - 1	2.6	30 - 1	3.64
31 - 1	2.56	31 - 1	3.7

Westgard	تفسیر	Creatinine	آزمایش	اتوآنالایزر BS-800	دستگاه	3	ماه	1402	سال
Normal Control	biorex	Mean	2.5	SD	0.1	Exp. Date	1403/1/17	Lot. No	035

