

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.11	
Mean	2.57	Mean	3.7
SD	0.05	SD	0.08
CV	1.8	CV	2.11
Day-Run	Value	Day-Run	Value
1 - 1	2.55	1 - 1	3.72
2 - 1	2.54	2 - 1	3.67
7 - 1	2.6	7 - 1	3.7
8 - 1	2.55	8 - 1	3.66
9 - 1	2.52	9 - 1	3.83
11 - 1	2.56	11 - 1	3.78
12 - 1	2.57	12 - 1	3.65
14 - 1	2.6	14 - 1	3.71
15 - 1	2.5	15 - 1	3.6
16 - 1	2.54	16 - 1	3.75
17 - 1	2.6	17 - 1	3.77
18 - 1	2.56	18 - 1	3.77
19 - 1	2.51	19 - 1	3.73
21 - 1	2.58	21 - 1	3.75
22 - 1	2.61	22 - 1	3.51
23 - 1	2.48	23 - 1	3.57
24 - 1	2.57	24 - 1	3.77
25 - 1	2.63	25 - 1	3.76
26 - 1	2.63	26 - 1	3.69
28 - 1	2.67	28 - 1	3.62
29 - 1	2.6	29 - 1	3.65
30 - 1	2.55	30 - 1	3.75

Westgard	تفسیر	Creatinine	آزمایش	BS-800 اتوآنالایزر	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	2.55	SD	0.1	Exp. Date	1403/1/17	Lot. No	035

