

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	0.45
SD	0.03
CV	7.36
Day-Run	Value
1 - 1	0.45
3 - 1	0.41
4 - 1	0.46
5 - 1	0.41
6 - 1	0.5
10 - 1	0.47
11 - 1	0.41
13 - 1	0.48
19 - 1	0.45
20 - 1	0.47
21 - 1	0.48
24 - 1	0.48
25 - 1	0.46
26 - 1	0.47
27 - 1	0.4
28 - 1	0.49
29 - 1	0.41
31 - 1	0.41

High - CV : 4.84	
Mean	2.2
SD	0.11
CV	4.84
Day-Run	Value
1 - 1	2.2
3 - 1	2.11
4 - 1	2.3
5 - 1	2.05
6 - 1	2.3
10 - 1	2.11
11 - 1	2.21
13 - 1	2.15
19 - 1	2.35
20 - 1	2.19
21 - 1	2.22
24 - 1	2.3
25 - 1	2.33
26 - 1	2.05
27 - 1	2.38
28 - 1	2.11
31 - 1	2.11

Westgard	تفسیر	D.Bili	آزمایش	BS-800 اتوآنالایزر	دستگاه	4	ماه	1402	سال
Normal Control	biorex	Mean	0.46	SD	0.03	Exp. Date	1403/1/17	Lot. No	035

