

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 5.8	
Mean	0.46	Mean	2.23
SD	0.03	SD	0.13
CV	6.15	CV	5.8
Day-Run	Value	Day-Run	Value
1 - 1	0.4	1 - 1	2.4
2 - 1	0.45	2 - 1	2.05
3 - 1	0.45	3 - 1	2.36
6 - 1	0.48	6 - 1	2.05
8 - 1	0.45	8 - 1	2.25
9 - 1	0.41	9 - 1	2.36
10 - 1	0.5	10 - 1	2
11 - 1	0.47	11 - 1	2.25
17 - 1	0.46	17 - 1	2.2
20 - 1	0.47	20 - 1	2.3
20 - 1	0.47	21 - 1	2.36
21 - 1	0.4	22 - 1	2.05
22 - 1	0.48	23 - 1	2.25
23 - 1	0.45	24 - 1	2.25
24 - 1	0.42	25 - 1	2.3
25 - 1	0.48	27 - 1	2.38
27 - 1	0.48	28 - 1	2.15
28 - 1	0.45	29 - 1	2.05
30 - 1	0.46	30 - 1	2.33
31 - 1	0.47	31 - 1	2.3

Westgard	تفسیر	D.Bili	آزمایش	BS-800 اتوآنالایزر	دستگاه	3	ماه	1402	سال
Normal Control	biorex	Mean	0.45	SD	0.03	Exp. Date	1403/1/17	Lot. No	035

