

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.85	
Mean	74.17	Mean	228.56
SD	1.25	SD	4.23
CV	1.68	CV	1.85
Day-Run	Value	Day-Run	Value
1 - 1	76	1 - 1	227
3 - 1	76	3 - 1	228
4 - 1	74	4 - 1	229
5 - 1	75	5 - 1	232
6 - 1	74	6 - 1	230
10 - 1	73	10 - 1	230
11 - 1	73	11 - 1	231
13 - 1	73	13 - 1	229
19 - 1	73	19 - 1	234
20 - 1	76	20 - 1	225
21 - 1	76	21 - 1	235
24 - 1	73	24 - 1	227
25 - 1	74	25 - 1	225
26 - 1	73	26 - 1	235
27 - 1	73	27 - 1	218
28 - 1	75	28 - 1	226
29 - 1	75	29 - 1	229
31 - 1	73	31 - 1	224

Westgard	تفسیر	FBS	آزمایش	BS-800 اتوآنالایزر	دستگاه	4	ماه	1402	سال
Normal Control	biorex	Mean	74	SD	1.39	Exp. Date	1403/1/17	Lot. No	035

