

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.48	
Mean	73.73	Mean	229
SD	0.83	SD	3.39
CV	1.12	CV	1.48
Day-Run	Value	Day-Run	Value
1 - 1	75	1 - 1	224
2 - 1	73	2 - 1	227
7 - 1	74	7 - 1	229
8 - 1	73	8 - 1	227
9 - 1	73	9 - 1	230
11 - 1	75	11 - 1	233
12 - 1	74	12 - 1	230
14 - 1	74	14 - 1	223
15 - 1	73	15 - 1	224
16 - 1	73	16 - 1	224
17 - 1	75	17 - 1	226
18 - 1	74	18 - 1	230
19 - 1	75	19 - 1	228
21 - 1	73	21 - 1	228
22 - 1	73	22 - 1	231
23 - 1	73	23 - 1	230
24 - 1	74	24 - 1	234
25 - 1	73	25 - 1	230
26 - 1	75	26 - 1	232
28 - 1	74	28 - 1	233
29 - 1	73	29 - 1	230
30 - 1	75	30 - 1	235

Westgard	تفسیر	FBS	آزمایش	BS-800 اتوآنالایزر	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	74.17	SD	1.25	Exp. Date	1403/1/17	Lot. No	035

