

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	22.64	
2 - 1	22.65	
3 - 1	22.69	
4 - 1	22.59	
7 - 1	22.72	
8 - 1	22.65	Cusum 3S
9 - 1	22.61	
10 - 1	22.65	
11 - 1	22.68	
12 - 1	22.69	
14 - 1	22.58	
15 - 1	22.62	
16 - 1	22.65	
17 - 1	22.54	
18 - 1	22.58	
19 - 1	22.65	
21 - 1	22.65	
22 - 1	22.67	
23 - 1	22.68	
24 - 1	22.64	
25 - 1	22.6	
26 - 1	22.62	
28 - 1	22.67	
29 - 1	22.61	
30 - 1	22.64	
31 - 1	22.67	

Westgard	تفسیر	FSH	آزمایش	2000 ایمولایت XPI-	دستگاه	5	ماه	1402	سال
Control	biorex	Mean	22.63	SD	0.05	Exp. Date	1403/4/2	Lot. No	365-02-028
Next Month	Mean	22.64	SD	0.04	CV	0.18			

