

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	7.52	
3 - 1	7.52	
4 - 1	7.65	
5 - 1	7.59	
6 - 1	7.45	
7 - 1	7.4	
10 - 1	7.34	
11 - 1	7.62	
12 - 1	7.59	
13 - 1	7.64	
14 - 1	7.45	
15 - 1	7.38	
17 - 1	7.61	
18 - 1	7.63	
19 - 1	7.45	
20 - 1	7.34	
21 - 1	7.45	
22 - 1	7.6	
24 - 1	7.52	
25 - 1	7.43	
26 - 1	7.39	
27 - 1	7.4	
28 - 1	7.46	
29 - 1	7.62	
31 - 1	7.6	

Westgard	تفسیر	Free T3	آزمایش	SIEMENS Centaur CP	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	7.48	SD	0.09	Exp. Date	1402/5/11	Lot. No	36-02-026
Next Month	Mean	7.51	SD	0.1	CV	1.36			

