

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

| Day-Run | Value | Description |
|---------|-------|-------------|
| 1 - 1 | 2.8 | |
| 2 - 1 | 2.94 | |
| 3 - 1 | 2.81 | |
| 4 - 1 | 2.65 | |
| 7 - 1 | 2.71 | |
| 8 - 1 | 2.85 | |
| 9 - 1 | 2.94 | |
| 10 - 1 | 2.59 | |
| 11 - 1 | 2.9 | |
| 12 - 1 | 2.75 | |
| 14 - 1 | 2.79 | |
| 15 - 1 | 2.95 | |
| 16 - 1 | 2.84 | |
| 17 - 1 | 2.75 | |
| 18 - 1 | 2.69 | |
| 19 - 1 | 2.78 | |
| 21 - 1 | 2.89 | |
| 22 - 1 | 2.84 | |
| 23 - 1 | 2.91 | |
| 24 - 1 | 2.8 | |
| 25 - 1 | 2.79 | |
| 26 - 1 | 2.75 | |
| 28 - 1 | 2.85 | |
| 29 - 1 | 2.92 | |
| 30 - 1 | 2.85 | |
| 31 - 1 | 2.88 | |

| | | | | | | | | | |
|------------|--------|---------|--------|--------------------|--------|-----------|------------|---------|-----------|
| Westgard | تفسیر | Free T4 | آزمایش | 2000 ایمولایت XPI- | دستگاه | 5 | ماه | 1402 | سال |
| Control | biorex | Mean | 2.81 | SD | 0.12 | Exp. Date | 1402/10/30 | Lot. No | 36-02-028 |
| Next Month | Mean | 2.82 | SD | 0.09 | CV | 3.27 | | | |

