

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	69	1:2S
3 - 1	73.3	
4 - 1	72.2	
5 - 1	72.1	
6 - 1	72.1	
10 - 1	70.9	
11 - 1	71.6	
13 - 1	70.5	
19 - 1	72.2	
20 - 1	69.6	
21 - 1	67.6	1:2S
24 - 1	73	
25 - 1	71.5	
26 - 1	70.4	
27 - 1	70.1	
28 - 1	68.5	
29 - 1	70.4	
31 - 1	70.7	

Westgard	تفسیر	Gamma GT	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	71.4	SD	1.45	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	70.87	SD	1.54	CV	2.18			

