

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	12.09
SD	0.21
CV	1.7
Day-Run	Value
1 - 1	12.1
3 - 1	12.5
4 - 1	12.2
5 - 1	11.8
6 - 1	12.1
7 - 1	11.8
10 - 1	12.1
11 - 1	12
12 - 1	11.9
13 - 1	12
14 - 1	12
15 - 1	12.1
16 - 1	12
17 - 1	12.4
18 - 1	12.3
19 - 1	12.3
20 - 1	12.2
21 - 1	11.7
22 - 1	12
24 - 1	12.2
25 - 1	11.9
26 - 1	12
27 - 1	12
28 - 1	12.5
29 - 1	12.2

High - CV : 2.04	
Mean	17.38
SD	0.35
CV	2.04
Day-Run	Value
1 - 1	17.3
3 - 1	17.4
4 - 1	17.2
5 - 1	17.2
6 - 1	17.1
7 - 1	17.2
10 - 1	17.4
11 - 1	17.1
12 - 1	17.1
13 - 1	17.1
14 - 1	17.2
15 - 1	18.3
16 - 1	18.2
17 - 1	17
18 - 1	18
19 - 1	17.3
21 - 1	17.1
22 - 1	17.2
24 - 1	17.2
25 - 1	17.5
26 - 1	17.3
27 - 1	17.4
28 - 1	17.8
29 - 1	17.6

Westgard	تفسیر	HB	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	4	ماه	1402	سال
Normal Control		Mean	12.2	SD	0.22	Exp. Date	1402/4/13	Lot. No	

