

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	39.72
SD	0.68
CV	1.72
Day-Run	Value
1 - 1	40
3 - 1	39.8
4 - 1	41.1
5 - 1	38.4
6 - 1	39.3
7 - 1	40
10 - 1	39.9
11 - 1	39.2
12 - 1	38.9
13 - 1	38.8
14 - 1	39.7
15 - 1	39.5
16 - 1	39.7
17 - 1	40
18 - 1	39.52
20 - 1	39.1
21 - 1	38.8
22 - 1	40.4
24 - 1	39.8
25 - 1	41
26 - 1	40.8
27 - 1	40
28 - 1	40
29 - 1	39.5

High - CV : 2.15	
Mean	54.82
SD	1.18
CV	2.15
Day-Run	Value
1 - 1	54.7
3 - 1	54.5
4 - 1	54.6
5 - 1	54.1
6 - 1	54.4
7 - 1	54.1
10 - 1	54.2
11 - 1	54.6
12 - 1	54.6
13 - 1	54.6
14 - 1	53.8
15 - 1	59
16 - 1	55
17 - 1	55.9
18 - 1	56
21 - 1	52.9
22 - 1	55.2
24 - 1	53.4
25 - 1	55.3
26 - 1	54.3
27 - 1	55.6
28 - 1	55.5
29 - 1	54.5

Westgard	تفسیر	HCT	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	4	ماه	1402	سال
Normal Control		Mean	39.58	SD	0.75	Exp. Date	1402/4/13	Lot. No	

