

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	39.76
SD	0.5
CV	1.27
Day-Run	Value
2 - 1	39
2 - 2	40
3 - 1	40
4 - 1	39
8 - 1	40.3
9 - 1	40
11 - 1	39.5
12 - 1	40
15 - 1	39.4
17 - 1	40
18 - 1	39.3
21 - 1	40.8
22 - 1	39.3
23 - 1	40.2
24 - 1	39
25 - 1	40.2
28 - 1	40
29 - 1	40
31 - 1	39.5

High - CV : 1.5	
Mean	54.94
SD	0.82
CV	1.5
Day-Run	Value
2 - 1	54.7
2 - 2	55
3 - 1	54.2
4 - 1	53.7
8 - 1	54.2
9 - 1	55.2
11 - 1	54.2
12 - 1	54.5
15 - 1	55.3
17 - 1	55.9
18 - 1	55.5
21 - 1	54.2
22 - 1	53.5
23 - 1	56
24 - 1	56.5
25 - 1	55
28 - 1	56
29 - 1	55.2
31 - 1	55

Westgard	تفسیر	HCT	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	5	ماه	1402	سال
Normal Control	Labex	Mean	39.72	SD	0.68	Exp. Date	1402/6/30	Lot. No	2305n

