

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	68	
3 - 1	66.1	
4 - 1	66.6	
5 - 1	66.5	
6 - 1	65	
10 - 1	66.7	
11 - 1	65.9	Cusum 3S
13 - 1	66.3	Cusum 3S
19 - 1	66.8	
20 - 1	67.6	
21 - 1	65.8	
24 - 1	70	1:2S
25 - 1	66	
26 - 1	67.6	
27 - 1	68	
28 - 1	65	
29 - 1	67	
31 - 1	64.6	

Westgard	تفسیر	HDL	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	67.08	SD	1.24	Exp. Date	1402/5/11	Lot. No	035
Next Month	Mean	66.64	SD	1.3	CV	1.95			

