

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	68.2	
2 - 1	65.9	
3 - 1	66.8	
6 - 1	66	
8 - 1	65.5	
9 - 1	65.3	
10 - 1	66.8	
11 - 1	67.5	
17 - 1	70	1:2S
20 - 1	67.1	
21 - 1	66.7	
22 - 1	68.1	
23 - 1	68.9	
24 - 1	67.6	
25 - 1	66.2	
27 - 1	65	
28 - 1	67	
29 - 1	67.2	
30 - 1	67.9	
31 - 1	67.8	

Westgard	تفسیر	HDL	آزمایش	اتوآنالایزر BS-800	دستگاه	3	ماه	1402	سال
Control	biorex	Mean	66	SD	1.45	Exp. Date	1402/10/2	Lot. No	035
Next Month	Mean	67.08	SD	1.24	CV	1.86			

