

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	127.5	
3 - 1	127.8	
4 - 1	122	
6 - 1	117.4	
7 - 1	120.4	
10 - 1	126.8	
11 - 1	127.6	
12 - 1	119.8	
13 - 1	119.6	
14 - 1	125.7	
15 - 1	128.6	
17 - 1	123.9	
18 - 1	121.5	
18 - 2	123	
20 - 1	125.4	
21 - 1	121.7	
22 - 1	127.6	
24 - 1	125.5	
25 - 1	121.6	
26 - 1	118.4	
28 - 1	129	
28 - 2	123.8	
29 - 1	122.3	
31 - 1	123.6	

Westgard	تفسیر	Ig A	آزمایش	BS-480 اتوآنالایزر	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	123.42	SD	3.2	Exp. Date	1402/5/11	Lot. No	641-01-030/2
Next Month	Mean	123.77	SD	3.38	CV	2.73			

