

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	125.9	
2 - 1	126.2	
3 - 1	124.6	
4 - 1	126.3	
6 - 1	126.2	
7 - 1	118.5	
8 - 1	127	
9 - 1	121.8	
10 - 1	122.3	
11 - 1	119	
13 - 1	120	
16 - 1	123.8	
17 - 1	125.9	
18 - 1	127	
20 - 1	118	
21 - 1	118.8	
22 - 1	127.2	
23 - 1	120.1	
24 - 1	129.5	
25 - 1	124.7	
27 - 1	123.6	
28 - 1	121.7	
29 - 2	121.3	
30 - 1	122.9	
31 - 1	123.1	

Westgard	تفسیر	Ig A	آزمایش	BS-480 اتوآنالایزر	دستگاه	3	ماه	1402	سال
Control	biorex	Mean	124.66	SD	3.64	Exp. Date	1402/1/9	Lot. No	641-01-030/2
Next Month	Mean	123.42	SD	3.2	CV	2.59			

