

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	557.6	
3 - 1	541.7	
3 - 2	538.2	
5 - 1	539.5	
6 - 1	537.2	
7 - 1	557.5	
10 - 1	523.8	
11 - 1	516.3	
12 - 1	533.8	
13 - 1	555.3	
14 - 1	526.8	
15 - 1	528	
17 - 1	557.5	
18 - 1	510.7	
20 - 1	539.6	
20 - 2	547.3	
21 - 1	537.2	
22 - 1	561.9	
24 - 1	508.8	
25 - 1	537.9	
26 - 1	528	
27 - 1	545	
29 - 1	536	
29 - 2	530.2	
31 - 1	564.4	

Westgard	تفسیر	Ig G	آزمایش	BS-480 اتوآنالایزر	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	537.6	SD	15.89	Exp. Date	1402/5/1	Lot. No	641-01-030/2
Next Month	Mean	538.41	SD	15.2	CV	2.82			

