

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	45.2	
3 - 1	46.4	
4 - 1	44.3	
6 - 1	45.6	
6 - 2	47.2	
7 - 1	47.7	
10 - 1	45.4	
11 - 1	48	
12 - 1	47.1	
13 - 1	48.2	
14 - 1	42.4	
15 - 1	44.2	
17 - 1	49.4	
18 - 1	41	
19 - 1	44.3	
20 - 1	46.3	
21 - 1	41.7	
22 - 1	45	
24 - 1	46	
25 - 1	44.9	
26 - 1	45.5	
27 - 1	48.4	
28 - 1	47.7	
29 - 1	43.2	
31 - 1	43.9	

Westgard	تفسیر	Ig M	آزمایش	BS-480 اتوآنالایزر	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	44.79	SD	2.42	Exp. Date	1402/5/11	Lot. No	641-01-030/2
Next Month	Mean	45.56	SD	2.15	CV	4.73			

