

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	237	
3 - 1	234	
4 - 1	228	
5 - 1	229.2	
6 - 1	233	
10 - 1	225.85	1:2S
11 - 1	225	1:2S 2:2S
13 - 1	228.5	
19 - 1	234.5	
20 - 1	225	1:2S
21 - 1	232	
24 - 1	229.5	
25 - 1	234	
26 - 1	236	
27 - 1	230.5	
28 - 1	233	
29 - 1	232.3	
31 - 1	233.1	

Westgard	تفسیر	Iron	آزمایش	BS-800 اتوآنالایزر	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	230	SD	3.14	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	231.14	SD	3.64	CV	1.58			

