

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	240	1:2S
2 - 1	232	
7 - 1	237.2	
8 - 1	230	
9 - 1	238.5	1:2S
11 - 1	232.6	
12 - 1	227	
14 - 1	233	
15 - 1	234	
16 - 1	236	
17 - 1	233.4	
18 - 1	238	
19 - 1	234	
21 - 1	233	
22 - 1	225	
23 - 1	223.2	1:2S
24 - 1	226.4	
25 - 1	235.1	4:1S
26 - 1	240	1:2S
28 - 1	229.5	
29 - 1	235	
30 - 1	232.1	

Westgard	تفسیر	Iron	آزمایش	BS-800 اتوآنالایزر	دستگاه	5	ماه	1402	سال
Control	biorex	Mean	231.14	SD	3.64	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	232.95	SD	4.65	CV	1.99			

