

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	137.5	
3 - 1	138	
4 - 1	133.9	
5 - 1	136.8	
6 - 1	137.8	
10 - 1	136.5	
11 - 1	137.2	
13 - 1	139	
19 - 1	140	
20 - 1	135	
21 - 1	133	1:2S
24 - 1	140	
25 - 1	140	
26 - 1	138	
27 - 1	132	1:2S
28 - 1	140	
29 - 1	134	
31 - 1	136	

Westgard	تفسیر	LDL	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	137.16	SD	2.5	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	136.93	SD	2.52	CV	1.84			

