

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	140	
2 - 1	133.9	
3 - 1	137.1	
6 - 1	135.4	
8 - 1	135	
10 - 1	138.8	
11 - 1	139.2	
17 - 1	136.1	
20 - 1	137.4	
21 - 1	137	
22 - 1	139.5	
23 - 1	139.9	
24 - 1	142	1:2S
25 - 1	138.9	
27 - 1	140	
28 - 1	140.9	
29 - 1	140	
30 - 1	136	
31 - 1	138	

Westgard	تفسیر	LDL	آزمایش	اتوآنالایزر BS-800	دستگاه	3	ماه	1402	سال
Control	biorex	Mean	137.53	SD	2.03	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	138.16	SD	2.2	CV	1.59			

