

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	42	
3 - 1	42.8	
4 - 1	43	
5 - 1	43.7	
6 - 1	42.5	
7 - 1	43	
10 - 1	44	
11 - 1	43.2	
12 - 1	42	
13 - 1	42.5	
14 - 1	43.2	
15 - 1	44	
17 - 1	42.5	
18 - 1	43.3	
19 - 1	43.8	
20 - 1	42.6	
21 - 1	43.5	
22 - 1	44	
24 - 1	42.6	
25 - 1	43.3	
26 - 1	44	
27 - 1	43.2	
28 - 1	42.5	
29 - 1	43	
29 - 2	44.1	
31 - 1	42	

Westgard	تفسیر	Lactate(Lactic Acid)	آزمایش	BS-480 اتوآنالایزر	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	43.13	SD	0.7	Exp. Date	1402/5/11	Lot. No	
Next Month	Mean	43.09	SD	0.66	CV	1.54			

