

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	43.1	
2 - 1	44	
3 - 1	42	
4 - 1	43.3	
6 - 1	42	
7 - 1	43.5	
8 - 1	43	
9 - 1	42.5	
10 - 1	43.5	
11 - 1	42.1	
13 - 1	42.9	
16 - 1	44	
17 - 1	43.1	
18 - 1	42.5	
20 - 1	43.6	
21 - 1	42	
22 - 1	44	
23 - 1	43.2	
24 - 1	43	
27 - 1	44.1	
28 - 1	42.5	
29 - 1	44	
30 - 1	43.1	
31 - 1	44	

Westgard	تفسیر	Lactate(Lactic Acid)	آزمایش	BS-480 اتوآنالایزر	دستگاه	3	ماه	1402	سال
Control	biorex	Mean	43.04	SD	0.82	Exp. Date	1401/11/18	Lot. No	
Next Month	Mean	43.13	SD	0.7	CV	1.63			

