

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	65.8	
2 - 1	68.9	
3 - 1	69.8	
6 - 1	68	
8 - 1	63.4	
9 - 1	63.2	
10 - 1	69	
11 - 1	69	
17 - 1	66	
20 - 1	66.5	
21 - 1	68.9	
22 - 1	70	
23 - 1	63.4	
24 - 1	68	
25 - 1	68	
27 - 1	68	
28 - 1	63.6	
29 - 1	66	
30 - 1	69.7	
31 - 1	63.7	

Westgard	تفسیر	Lipase	آزمایش	اتوآنالایزر BS-800	دستگاه	3	ماه	1402	سال
Control	biorex	Mean	65.3	SD	2.51	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	66.95	SD	2.4	CV	3.59			

