

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	69	
2 - 1	65	
7 - 1	63.7	
8 - 1	69.7	
9 - 1	68	
11 - 1	65	
12 - 1	63.8	
14 - 1	69.4	
15 - 1	65.1	
16 - 1	66	
17 - 1	69.7	
18 - 1	69.1	
19 - 1	68.2	
21 - 1	69	
22 - 1	67.5	
23 - 1	64.8	
24 - 1	69	
25 - 1	67.2	
26 - 1	68	
28 - 1	67.5	
29 - 1	68	
30 - 1	69	

Westgard	تفسیر	Lipase	آزمایش	BS-800 اتوآنالایزر	دستگاه	5	ماه	1402	سال
Control	biorex	Mean	67.25	SD	2.36	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	67.35	SD	1.97	CV	2.93			

