

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

| Normal | |
|---------|-------|
| Mean | 28.16 |
| SD | 0.44 |
| CV | 1.55 |
| Day-Run | Value |
| 1 - 1 | 27.9 |
| 2 - 1 | 28.2 |
| 4 - 1 | 27.5 |
| 5 - 1 | 28 |
| 6 - 1 | 28 |
| 7 - 1 | 29 |
| 10 - 1 | 28.3 |
| 11 - 1 | 28.3 |
| 12 - 1 | 28.7 |
| 13 - 1 | 28 |
| 14 - 1 | 27.7 |
| 15 - 1 | 29.2 |
| 17 - 1 | 28 |
| 18 - 1 | 28.8 |
| 19 - 1 | 28 |
| 20 - 1 | 28.2 |
| 18 - 2 | 28.27 |
| 20 - 2 | 28.3 |
| 21 - 1 | 27.6 |
| 22 - 1 | 28 |
| 24 - 1 | 28.8 |
| 25 - 1 | 27.8 |
| 26 - 1 | 27.8 |
| 27 - 1 | 28.1 |
| 28 - 1 | 28.1 |
| 29 - 1 | 27.6 |

| High - CV : 1.05 | |
|------------------|-------|
| Mean | 31.56 |
| SD | 0.33 |
| CV | 1.05 |
| Day-Run | Value |
| 1 - 1 | 31.4 |
| 2 - 1 | 31.7 |
| 4 - 1 | 31.7 |
| 5 - 1 | 31.9 |
| 6 - 1 | 31.6 |
| 7 - 1 | 31.6 |
| 10 - 1 | 32 |
| 11 - 1 | 31.6 |
| 12 - 1 | 31.4 |
| 13 - 1 | 31.2 |
| 14 - 1 | 31.2 |
| 15 - 1 | 31.5 |
| 17 - 1 | 31 |
| 18 - 1 | 31.3 |
| 19 - 1 | 32 |
| 20 - 1 | 31.7 |
| 18 - 2 | 31.2 |
| 21 - 1 | 31.6 |
| 22 - 1 | 30.9 |
| 24 - 1 | 32.2 |
| 25 - 1 | 31.8 |
| 26 - 1 | 31.9 |
| 27 - 1 | 31.6 |
| 28 - 1 | 31.9 |
| 29 - 1 | 31.2 |

| | | | | | | | | | |
|----------------|-------|------|--------|-----------------------|--------|-----------|-----------|---------|-----|
| Westgard | تفسیر | MCH | آزمایش | NIHON KOHDEN سل کانتر | دستگاه | 4 | ماه | 1402 | سال |
| Normal Control | | Mean | 28.26 | SD | 0.47 | Exp. Date | 1402/4/13 | Lot. No | |

