

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	30.6	
2 - 1	30.2	
3 - 1	30.9	
4 - 1	30	
5 - 1	31.2	
7 - 1	30	
8 - 1	31.3	
9 - 1	31.2	
10 - 1	30.6	
11 - 1	29.5	
12 - 1	29.8	
16 - 1	30.6	
17 - 1	30.5	
18 - 1	30	
20 - 1	30.9	
21 - 1	31.1	
23 - 1	30.8	
24 - 1	31.5	
25 - 1	30.6	
27 - 1	31	

Westgard	تفسیر	MCHC	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	3	ماه	1402	سال
Control	Labex	Mean	30.56	SD	0.55	Exp. Date	1402/2/9	Lot. No	
Next Month	Mean	30.62	SD	0.55	CV	1.8			

